

b

Double-Lumen Endobronchial Tube

Made from non -toxic PVC, 100% latex free

Available with left-sided and right-sided types

The cuffs, two lumens and pilot balloons are imprinted with tracheal (clear) or bronchial (blue) for easy identification

Smooth inner lumens allow easy passage of bronchoscopes

Pre-loaded with stylet to help tube maintain shape for easier intubation

Packed with switch connector and two low friction suction catheters



Left Side

REF No.	Size (Fr)
131028	28
131032	32
131035	35
131037	37
131039	39
131041	41

Right Side

REF No.	Size (Fr)
132028	28
132032	32
132035	35
132037	37
132039	39
132041	41