

Endotracheal Tube Standard high volume low pressure cuff

Made from non-toxic PVC, kink-resistant tube

100% latex free

Softer rounded Murphy eye is less invasive

Accurate depth marks

Standard connector 15mm

Radio opaque line provided



REF No.	Size ID (mm)
101040	4.0
101045	4.5
101050	5.0
101055	5.5
101060	6.0
101065	6.5
101070	7.0
101075	7.5
101080	8.0
101085	8.5
101090	9.0
101095	9.5
101100	10.0

Endotracheal Tube Standard low profile cuff

Made from non-toxic PVC, kink-resistant tube

100% latex free

Low profile cuff for short-term surgical intubation

Softer rounded Murphy eye is less invasive

Accurate depth marks

Standard connector 15mm

Radio opaque line provided



REF No.	Size ID (mm)
102040	4.0
102045	4.5
102050	5.0
102055	5.5
102060	6.0
102065	6.5
102070	7.0
102075	7.5
102080	8.0
102085	8.5
102090	9.0
102095	9.5
102100	10.0

Endotracheal Tube without Cuff

Made from non-toxic PVC, kink-resistant tube

100% latex free

Softer rounded Murphy eye is less invasive

Accurate depth marks

Standard connector 15mm

Radio opaque line provided



REF No.	Size ID (mm)
103030	3.0
103035	3.5
103040	4.0
103045	4.5
103050	5.0
103055	5.5
103060	6.0
103065	6.5
103070	7.0
103075	7.5
103080	8.0
103085	8.5
103090	9.0
103095	9.5
103100	10.0

Oral Preformed Tube High Volume low pressure cuff

Made from non-toxic PVC, kink-resistant tube

100% latex free

Softer rounded Murphy eye is less invasive

Accurate depth marks

Standard connector 15mm

Radio opaque line provided



REF No.	Size ID (mm)
104040	4.0
104045	4.5
104050	5.0
104055	5.5
104060	6.0
104065	6.5
104070	7.0
104075	7.5
104080	8.0
104085	8.5
104090	9.0
104095	9.5
104100	10.0

Oral Preformed Tube low profile cuff

Made from non-toxic PVC, kink-resistant tube

100% latex free

Low profile cuff for short-term surgical intubation

Softer rounded Murphy eye is less invasive

Accurate depth marks

Standard connector 15mm

Radio opaque line provided



REF No.	Size ID (mm)
105040	4.0
105045	4.5
105050	5.0

105055	5.5
105060	6.0
105065	6.5
105070	7.0
105075	7.5
105080	8.0
105085	8.5
105090	9.0
105095	9.5
105100	10.0

Oral Preformed Tube uncuffed

Made from non-toxic PVC, kink-resistant tube

100% latex free

Softer rounded Murphy eye is less invasive

Accurate depth marks

Standard connector 15mm

Radio opaque line provided



REF No.	Size ID (mm)
106030	3.0
106035	3.5
106040	4.0
106045	4.5
106050	5.0
106055	5.5
106060	6.0
106065	6.5
106070	7.0
106075	7.5
106080	8.0
106085	8.5
106090	9.0
106095	9.5
106100	10.0

Nasal Preformed Tube High Volume low pressure cuff

Made from non-toxic PVC, kink-resistant tube

100% latex free

Softer rounded Murphy eye is less invasive

Accurate depth marks

Standard connector 15mm

Radio opaque line provided



REF No.	Size ID (mm)
107040	4.0
107045	4.5
107050	5.0
107055	5.5
107060	6.0
107065	6.5
107070	7.0
107075	7.5
107080	8.0
107085	8.5
107090	9.0

Nasal Preformed Tube low profile cuff

Made from non-toxic PVC, kink-resistant tube

100% latex free

Low profile cuff for short-term surgical intubation

Softer rounded Murphy eye is less invasive

Accurate depth marks

Standard connector 15mm

Radio opaque line provided



REF No.	Size ID (mm)
108040	4.0
108045	4.5
108050	5.0
108055	5.5
108060	6.0
108065	6.5
108070	7.0
108075	7.5
108080	8.0
108085	8.5
108090	9.0
108095	9.5

108100

10.0

Nasal Preformed Tube uncuffed

Made from non-toxic PVC, kink-resistant tube

100% latex free

Softer rounded Murphy eye is less invasive

Accurate depth marks

Standard connector 15mm

Radio opaque line provided



REF No.	Size ID (mm)
109030	3.0
109035	3.5
109040	4.0
109045	4.5
109050	5.0
109055	5.5
109060	6.0
109065	6.5
109070	7.0
109075	7.5
109080	8.0
109085	8.5
109090	9.0
109095	9.5
109100	10.0

Reinforced Endotracheal Tube with cuff

Made from non-toxic PVC, kink-resistant tube

100% latex free

Flex conform to any patient positions, especially to OPS of decubitus

Softer rounded Murphy eye is less invasive

Accurate depth marks

Standard connector 15mm

Radio opaque line provided



REF No.	Size ID (mm)
141040	4.0
141045	4.5
141050	5.0
141055	5.5
141060	6.0
141065	6.5
141070	7.0
141075	7.5
141080	8.0
141085	8.5
141090	9.0

Reinforced Endotracheal Tube without cuff

Made from non-toxic PVC, kink-resistant tube

100% latex free

Flex conform to any patient positions, especially to OPS of decubitus

Softer rounded Murphy eye is less invasive

Accurate depth marks

Standard connector 15mm

Radio opaque line provided



REF No.	Size ID (mm)
142030	3.0
142035	3.5
142040	4.0
142045	4.5
142050	5.0
142055	5.5
142060	6.0
142065	6.5
142070	7.0
142075	7.5
142080	8.0
142085	8.5
142090	9.0
142095	9.5
142100	10.0