

# Personal Protective Equipment

## Face mask, 3ply, ear loops, TYPE II



3-layer material including:

PP Spunbonded + PP Meltblown + PP Spunbonded

BEF $\geq$ 98%

Differential Pressure(Pa/cm<sup>2</sup>): < 29.49 (Breathing Resistance)

Standard(s) applied:

EN 14683:2019, TYPE II

REF	Type	Color
2020001	Type II	Blue

## Face mask, 4ply, ear loops, TYPE IIR



99F1 - Type IIR

4-layer material including:

PP Spunbonded + PET Spunlaced + PP Meltblown + PP Spunbonded

BEF $\geq$ 99%


Differential Pressure(Pa/cm<sup>2</sup>): < 60 (Breathing Resistance)

Fluid resistance pressure (mmHg): > 120mmHg

Standard(s) applied: EN 14683:2019, TYPE IIR

REF	Type	Color
2020006	Type IIR	Blue

### Kids Disposable Face Mask 3-ply, ear loops



3-layer material including

PP Spunbonded + PP Meltblown + PP Spunbonded

BEF $\geq$ 98%

Differential Pressure(Pa/cm<sup>2</sup>): <29.49 (Breathing Resistance)

Standard(s) applied:

EN 14683:2019, TYPE II

REF	Type	Color
2020001-S	Type II	Carton print (white, blue, pink)