









Personal Protective Equipment

| | |
|--|---|
|  |  |
|  |  |
| <p align="center">Protective Goggles REF: LB008C</p> | <p align="center">Safety Glasses REF: LB006</p> |
| <p>Color: Clear</p> | <p>Color: Clear</p> |
| <p>76 g</p> | <p>39.6 g</p> |
| <p>Frame Material: PC</p> | <p>Frame Material: PC</p> |
| <p>Lens Material: Acetate</p> | <p>Lens Material: PC</p> |
| <p>With black elastic band and four vent holes with white sealing plugs</p> | |
| <p>Feature: anti-fog, anti-spittle</p> | <p>Feature: anti-fog, anti-spittle</p> |
| <p>Packing : 1pc/OPP bag, 15pcs/inner box 150pcs/master carton</p> | <p>Packing : 1pc/OPP bag, 20pcs/inner box 200pcs/master carton</p> |
| <p>Carton size: 41*41*68cm, 13.5kg</p> | <p>Carton size: 41*41*68cm, 10kg</p> |
| <p>FDA</p> | <p>FDA</p> |

| | |
|---|---|
|  |  |
| <p align="center">Safety Protective Goggles REF: LB008C</p> | <p align="center">Safety Protective Goggles REF: MT009</p> |
| <p>Color: Clear</p> | <p>Color: Clear</p> |
| <p>73.6 g</p> | <p>106 g</p> |
| <p>Frame Material: PC</p> | <p>Frame Material: TPU</p> |
| <p>Lens Material: Acetate</p> | <p>Lens Material: PC</p> |
| <p>With black elastic band</p> | <p>With black elastic band</p> |
| <p>Feature: anti-fog, anti-spittle</p> | <p>Feature: anti-fog, anti-spittle</p> |
| <p>Packing : 1pc/OPP bag, 15pcs/inner box 150pcs/master carton</p> | <p>Packing : 1pc/OPP bag, 15pcs/inner box 150pcs/master carton</p> |
| <p>Carton size: 52*40*51cm, 13.5kg</p> | <p>Carton size: 59.5*44.5*57cm, 17kg</p> |
| <p>FDA</p> | <p>FDA</p> |
|  |  |
| <p align="center">Safety Protective Goggles REF: JY001</p> | <p align="center">Protective visors face shield REF: LB009</p> |
| <p>Color: Clear</p> | <p>Color: Clear</p> |
| <p>75 g</p> | <p>40 g</p> |
| <p>Frame Material: TPU</p> | <p>Material: PET</p> |
| <p>Lens Material: PC</p> | <p>shield size: 33*22*0.025cm</p> |
| <p>With black elastic band</p> | <p>With white sponge and elastic strap</p> |
| <p>Feature: anti-fog, anti-spittle</p> | <p>Feature: Anti splash, anti-spittle</p> |
| <p>Packing: 1pc/OPP bag, 15pcs/inner box 150pcs/master carton</p> | <p>Packing: 5pc/OPP bag, 120pcs/ carton</p> |
| <p>Carton size: 41*41*68cm, 12.25kg</p> | <p>Carton size: 60*26*30cm, 8kg</p> |
| <p>FDA</p> | <p>FDA</p> |