


## Silicone Round Drainage



REF No.	Tube Size (FR)	Tube length (cm)	OD (mm)
124009	9	50	3.0
125012	12	50	4.0
125014	14	50	4.7
125016	16	50	5.3
125018	18	50	6.0
125020	20	50	6.7
125024	24	50	8.0
125026	26	50	8.7
125030	30	50	10.0
125033	33	50	11.0
125036	36	50	12.0

## Silicone Round Drainage Systems



- Round-drain 100 cm, made of 100% medical grade silicone , central opening with 4 lateral atraumatic eyes
- provides gentle drainage with gravity
- Continuous x-ray opaque line
- Available Bag 350 ml / 500 ml, with NRV and drain valve, incl. 2 holders
- Latex-free
- DEHP-free
- Shelf life: 5 years

REF No.	Tube Size (FR)	Bag Volume (ml)	Tube length (cm)	OD (mm)
114012	12	350	100	4.0
114014	14	350	100	4.7
114016	16	350	100	5.3
114018	18	350	100	6.0
114020	20	350	100	6.7
114024	24	350	100	8.0
114026	26	350	100	8.7
114030	30	350	100	10.0
114033	33	350	100	11.0
114036	36	350	100	12.0

REF No.	Tube Size (FR)	Bag Volume (ml)	Tube length (cm)	OD (mm)
115012	12	500	100	4.0
115014	14	500	100	4.7
115016	16	500	100	5.3
115018	18	500	100	6.0

115020	20	500	100	6.7
115024	24	500	100	8.0
115026	26	500	100	8.7
115030	30	500	100	10.0
115033	33	500	100	11.0
115036	36	500	100	12.0

## Silicone Soft Capillary Drainage



- Continuous x-ray opaque line
- Inner ribbing enables gentle capillary wound drainage and helps keep drain from clogging
- Latex-free
- Shelf life: 5 years

REF No.	Tube Size (FR)	Tube length (cm)	Diameter mm (mm)
121006	6	30	6
121008	8	30	8
121010	10	30	10
121012	12	30	12
121106	6	50	6
121108	8	50	8
121110	10	50	10
121112	12	50	12