

Urology

Urine Bag

2000ml Luxury Urine Bag



- 2,000 ML with air-vent
- Reinforced double hanger and cord with small hook,
- Sampling Port with Cap, and Big T-valve
- 90~120cm Inlet Tube with Bed Sheet Clip and big clamp
- Three-part Anti-Reflux Drip Chamber
- Silicone Outlet Tube with Tube Pocket
- Latex free
- Shelf life: 5 years

REF No.	Bag Volume (ml)	Tube length (cm)	Sampling Port	Valve
030090	2000	90	Stepping Sampling Port	Big T valve, horizontal
030090-N	2000	90	MONTO Needle free	Big T valve, horizontal
030092	2000	90	Stepping Sampling Port	Big T valve, vertical
030092-N	2000	90	MONTO Needle free	Big T valve, vertical
030210	2000	120	Stepping Sampling Port	Big T valve, horizontal
030210-N	2000	120	MONTO Needle free	Big T valve, horizontal
030212	2000	120	Stepping Sampling Port	Big T valve, vertical
030212-N	2000	120	MONTO Needle free	Big T valve, vertical

5000ml Luxury Urine Bag



- 5,000 ML with air-vent
- Reinforced double hanger and cord with small hook
- Sampling Port with Cap, and Big T-valve
- 90~120cm Inlet Tube with Bed Sheet Clip and big clamp
- Drip Chamber with Anti-Reflux
- Silicone Outlet Tube with Tube Pocket
- Latex free
- Shelf life: 5 years

REF No.	Bag Volume (ml)	Tube length (cm)	Sampling Port	Valve
030250	5000	120	Stepping Sampling Port	Big T valve

4000ml Luxury Urine Bag



- 4,000 ML with air-vent
- Reinforced double hanger and cord with small hook
- Sampling Port with Cap, and Big T-valve
- 90~120cm Inlet Tube with Bed Sheet Clip and big clamp
- Drip Chamber with Anti-Reflux
- Silicone Outlet Tube with Tube Pocket
- Latex free
- Shelf life: 5 years

REF No.	Bag Volume (ml)	Tube length (cm)	Sampling Port	Valve
030240	4000	120	MONTO Needle free	Big T valve

Luxury Urine Bag with Big Holder



- Luxury urine bag
- 2000 ML with air-vent
- With PVC non-return valve (NRV)
- Stepping Sampling Port and with Cap
- 90~120cm Inlet Tube with Bed Sheet Clip and big clamp
- With big holder
- Latex Free
- Shelf life: 5 years

REF No.	Bag Volume (ml)	Tube length (cm)	Sampling Port	Valve
030120-H	2000	120	Stepping Sampling Port	Big T valve

Luxury Urine Bag without Hanger



- 2,000 ML with Air-vent
- with PVC non-return valve(NRV)
- Stepping Sampling Port and with Cap
- 90~120cm Inlet Tube with Bed Sheet Clip and big clamp
- Latex Free
- Shelf life: 5 years

REF No.	Bag Volume (ml)	Tube length (cm)	Sampling Port	Valve
030190	2000	90	Stepping Sampling Port	Big T valve
030120	2000	120	Stepping Sampling Port	Big T valve

Economy Urine Bag without T valve



- With universal step connector
- Double sealed film
- Extra reinforced suspension holes
- Sterile version (EO) with non-return valve (NRV)
- Latex-free
- Shelf life: 5 years

REF No.	Bag Volume (ml)	Tube length (cm)	Valve
030615	1500	90	NA
030620	2000	90	NA

Economy Urine Bag with T valve



- With universal step connector
- Double sealed film
- Extra reinforced suspension holes
- Big T valve for easy one-hand drainage
- Non-return valve (NRV)
- Latex-free
- Shelf life: 5 years

REF No.	Bag Volume (ml)	Tube length (cm)	Inlet connector	Valve
030220	2000	90	Universal Step Connector	Big T valve

Biliary Drainage Bag



- With or without sterile,
- Without outlet and drainage valve
- With a special connector connected to different drainage system which is commonly used to drainage secretion and biliary.
- Made of PVC and latex-free rubber.
- Bag volume is 1,500ml, which can be customized

REF No.	Bag Volume (ml)	Tube length (cm)	Valve
030625	1500	90	NA
030630	1500	90	NA